

# CheckLite™ 250 Plus



Kikkoman Biochemifa Company

## Microbial Biomass Assay Kit



### Kit Features

- CheckLite™ 250 Plus is a microbial biomass assay kit based on the ATP measurement method.
- CheckLite™ contains thermostable firefly luciferase.
- CheckLite™ is characterized by its ease of use and fast performance.
- CheckLite™ can detect 10,000 coliform cells/mL (1,000 cells/assay).

### Kit Composition

<input type="checkbox"/> Luciferin-luciferase reagent (lyophilized)	5
<input type="checkbox"/> Reconstitution Buffer (5.5 ml)	5
<input type="checkbox"/> ATP-Releasing Reagent (5.5 mL)	5

### 250 Assay / Kit

CheckLite™ 250 Plus should be stored at a temperature between 2 and 8 °C in the refrigerator. **DO NOT FREEZE IT.**



### +CheckLite™ ATP Removal Kit

## HyServe

HyServe GmbH  
Ammerthalstraße 7  
Tel. +49 (0) 89 3703 1223  
info@hyserve.com

D - 85551  
Kirchheim/München  
Fax +49 (0) 89 3703 1225  
www.hyserve.com

# CheckLite™ HS Set

## High Sensitivity Microbial Biomass Assay Kit

**Luciferin-Luciferase  
ATP Assay Reagent  
Kit for  
Hypersensitive  
Viable Count Assay**



### Kit Features


- CheckLite™ HS Set is a microbial biomass assay kit based on the ATP measurement method.
- CheckLite™ contains thermostable firefly luciferase and ATP-scavenging enzymes.
- CheckLite™ is characterized by ease of use, fast performance, and high sensitivity.
- CheckLite™ can detect 1,000 coliform cells/mL (100 cells/assay).

### Kit Composition

- ☐ HS luciferin-luciferase reagent (lyophilized)
- ☐ Reconstitution Buffer (5.5 ml)
- ☐ ATP Scavenging Reagent (Freeze-Dried)
- ☐ Reconstitution Buffer (5.5 ml)
- ☐ ATP-Releasing Reagent (5.5 mL)

100 Assays/Kit

### Sample Essays Biomass



**The CheckLite™ HS kit should be stored  
at a temperature between 2 and 8 °C in  
the refrigerator. DO NOT FREEZE IT.**

**HyServe**

HyServe GmbH  
Ammerthalstraße 7  
Tel. +49 (0) 89 3703 1223  
[info@hyserve.com](mailto:info@hyserve.com)

D- 85551  
Kirchheim/München  
Fax +49 (0) 89 3703 1225  
[www.hyserve.com](http://www.hyserve.com)

**Kit for microbial biomass testing of aseptic products**

# CheckLite™ AT 100

CheckLite AT 100 is a suitable kit for time-saving microbial biomass testing on aseptic products, such as milk, cream, whipped cream, canned coffee, etc.



## Kit Features

- CheckLite® AT100 is a kit for highly sensitive microbial biomass assays based on the ATP measurement method.
- Kikkoman's ATP Removal Reagent and Treatment Reagent allow the centrifugation step, which is often required for conventional ATP methods, to be skipped.
- CheckLite® AT100 significantly reduces time and labor for aseptic testing of products such as milk, cream, whipped cream, canned coffee, etc.

## Kit Composition

<input type="checkbox"/> HS luciferin-luciferase reagent (lyophilized)	2
<input type="checkbox"/> Reconstitution Buffer for HS Luciferin-Luciferase Reagent (5.5 mL)	2
<input type="checkbox"/> ATP Scavenging Reagent (Freeze-Dried)	2
<input type="checkbox"/> Reconstitution buffer for ATP reagent	2
<input type="checkbox"/> ATP-Releasing Reagent	2
<input type="checkbox"/> Sample dilution buffer	4
<input type="checkbox"/> Sample Treatment Reagent	2

## 100 Essays/Kit

**CheckLite AT 100**  
should be stored at a  
temperature between  
**2 and 8 °C in the  
refrigerator. DO NOT  
FREEZE IT.**


**HyServe**

HyServe GmbH  
Ammerthalstraße 7  
Tel. +49 (0) 89 3703 1223  
info@hyserve.com

D- 85551  
Kirchheim/München  
Fax +49 (0) 89 3703 1225  
www.hyserve.com

## Procedure

### 1 Sample Preparation




Samples that have been incubated for a certain amount of time for an aseptic test can be used directly or diluted with buffer.

### 2 Sample Pretreatment



Place the sample treatment reagent in a test tube and mix it with the sample.

### 3 Non-microbial ATP removal



Add the sample dilution buffer and ATP remover solution into the test tube.

### 4 Microbial ATP extraction






Transfer the aliquot to another test tube and add the ATP-releasing reagent.

### 5 Luminescence measurement



Add the Luciferin-Luciferase HS reagent and measure the luminescence.

## Example



HyServe

HyServe GmbH  
Ammerthalstraße 7  
Tel. +49 (0) 89 3703 1223  
info@hyserve.com

D- 85551  
Kirchheim/München  
Fax +49 (0) 89 3703 1225  
www.hyserve.com