

Trypto-Soya Agar (SCD Agar)

(Soybean-Casein-Digest Agar)



Code	Form	Package	Storage	Shelf Life
05516	Granule	300 g	Dry, RT	3 years

Directions

Suspend 40.0 g of the dehydrated medium in 1,000 mL of distilled water, mix well and heat to dissolve the medium. Sterilize by autoclaving at 121 °C for 15 minutes. Use as plates or slants.

Remarks

The medium is fit for various purposes because it permits the growth of fastidious bacteria that do grow on Nutrient Agar or Heart Infusion Agar, and is markedly superior to conventional media as a basal medium of blood agar.

Formula

Components	In 40.0 g/L
Peptone	15.0 g
Soya Peptone	5.0 g
Sodium Chloride	5.0 g
Agar	15.0 g

pH 7.3 ± 0.2